

PRODUCT DATA SHEET for 2000-0785

4.5 mL Internally threaded and capped (includes o-ring)
Cryogenic vial, Stargrip base, Sterile

Product Specifications:

Packaging: 100 per bag, 10 bags per case

Sterility: Sterile (Up to 5 years)

Temperature:

- Suitable for a wide range of temperatures—from -196° to +121°C.
 - Optimized for vapor-phase liquid nitrogen storage.
- Note: not suitable for immersion in liquid nitrogen.

Dimensions: Height: 90.9mm; Diameter: 12.7mm (at widest point), 12.7mm cap diameter, 11.7mm internal diameter

Racking: Not Racked

Caps are assembled onto tubes.

Cap Specifications:

Color: Clear

Material: Polypropylene cap has silicone rubber O-ring for a secure, leak proof seal.

Cap Thread: Internally-Threaded

Seal: Silicone O-Ring

Tube Specifications:


Volume: 4.5 mL

Base: Stargrip

Color: Clear

Material:

- Ultra-pure polypropylene meets FDA 21 CFR 177.1520 & USP Class VI.
- High-contrast black 100 µL graduations are bold and easy to read.
- Large white writing area on all tubes.
- Maximize tube storage within racks, boxes, and canes.
- Sterile and certified free of detectable DNase, RNase, DNA, PCR inhibitors, and endotoxins.



The specifications provided herein are based on the information supplied by and tests conducted by primary material manufacturer (s). The end-user shall independently determine the suitability of the tubes for the purpose for which they were purchased.